


7 Use OMERO.web to view your data from any web browser with the URL:

nightshade.openmicroscopy.org

You do not have to enter a server name
 OMERO.web uses the same login, and has the same functionality, as OMERO.insight, except OMERO.web cannot import data



Detailed User Guides for common workflows using OMERO.insight and OMERO.web are available online and for download at:

help.openmicroscopy.org

For help to get started using OMERO or any other questions please email us at:

ux@openmicroscopy.org.uk

Getting Started with OMERO.insight and Nightshade

CLS has an OMERO server called Nightshade that allows you to store image data centrally. You can use the OMERO.insight client to upload, view and download data from any computer or laptop on or off-campus.

Your University of Dundee username and password are used to log into the OMERO server, using the URLs given below. Your OMERO account is automatically created the first time you log in and if you change your University of Dundee password, OMERO will automatically reflect this update.

1 Download the OMERO.insight client for your operating system from:

http://downloads.openmicroscopy.org/latest/omero5



OMERO 5.0.0 Downloads

[Clients](#) | [Plugins](#) | [Additional](#) | [Servers](#) | [Virtual Appliance](#) | [API](#) | [Code](#) | [Components](#) | [Legacy](#)


- Information on this release of OMERO is in the [release announcement](#)
- Full documentation is available as [web documentation](#) or [PDF documentation](#) and there are user guides for the clients on our [Help website](#)
- A standard OMERO user just needs to download the client package with the same major version as their institutional server e.g. 5.0 clients with the 5.0 server
- If you do not have an institutional server, you can apply for an [account on our Demo server](#) or download the [Virtual Appliance](#) to install your own version locally.

OMERO client downloads

| Clients | Size | File Name | Checksum |
|----------|----------|---|--------------|
| Windows | 70.2 MB | OMERO.clients-5.0.0-ice35-b19.win.zip | 002ce5 (MD5) |
| Mac OS X | 205.4 MB | OMERO.clients-5.0.0-ice35-b19.mac.zip | 625a57 (MD5) |
| Linux | 68.34 MB | OMERO.clients-5.0.0-ice35-b19.linux.zip | fe8426 (MD5) |

- Each client package includes [OMERO.insight](#), [OMERO.importer](#) and [OMERO.editor](#) and requires Java Version 1.6 or higher. OMERO.web is part of the server package, so individual users do not need to install it locally.
- Full instructions for installing the clients are on the Help website: [Getting Started with OMERO.insight Version 5.0.0](#)


2 Expand the OMERO.insight client archive




- 3 Open the OMERO.insight client
Click the **Spanner** icon



- 5 Enter your University of Dundee **Username** and **Password**
Click **Login**



- 4 In the **Add Server** dialogue box, click the **+** icon
Enter the server address:
nightshade.openmicroscopy.org.uk
Click **Apply**



- 6 Click the **Importer** icon to add image data
Using the file chooser in the Importer window select and queue data for import
Use the location window to choose existing or create new Projects, Datasets or Screens into which to import the data



Note: To import data from a networked drive on Macs, the drive must be mounted in the Finder first