

Using OMERO.insight and OMERO.web

OMERO stores image data on a central server. You can use the OMERO.insight client to upload, view and download data from any personal computer.

Your institution's OMERO server administrator will provide you with the server address to use when connecting from OMERO.insight. They may have set up a user name and password for you too. Alternatively your username and password may be your standard login for the institution's networked computers.

Installing

1 Download OMERO.insight Client

https://www.openmicroscopy.org/site/support/omero4/downloads



2 Download OMERO.insight Client







Open OMERO.insight Client

 OMERO.insight OMERO.insight.app 	Starting data management services
MEROinsight_unix.sh	howe.openmicroscopy.org.uk [LAN] 🔟 🖏
	Username: polly
	Password:
	Login
	4.40-0003-0003 OMERO is distributed under the terms of the GNU GPL. For more information, visit openmicroscopy.org
	••• OME





Adding Server and Logging In

Adding the server is only required the first time an OMERO server is accessed from the computer.

Click on th	e Spanner	icon.
-------------	------------------	-------

1

••••••••••••••••••••••••••••••••••••••	
localhost [LAN] 🛋 🖏	Z
Username:	
Password:	
Quit	
4.4.4252-787690-01th-1683-553 OMERO is distributed under the terms of the GNU GPL. For more information, visit openmicroscopy org	

2 In the Add Server dialogue box, click on the + icon.

	Servers	×
	Servers Enter a new server or select an existing one.	Z
	Server Address and Port e.g. test.openmicroscopy.org or 134.20.12.33 You should not have to modify the port.	
	Iocalhost	4064
(+)	Connection Speed	
	Apply	Cancel

••• OME



3 Enter the server address provided by your OMERO administrator. Click on **Apply**.



4 Enter the **Username** and **Password** provided by your OMERO administrator. Click on **Login**.

